



DO MORE WITH

FORA[®] 6

Connect



Stay informed about your health and make better wellness decisions with FORA[®] 6 Connect.



Blood Glucose
Hematocrit
Hemoglobin



β -Ketone



Uric Acid



Total
Cholesterol



REVOLUTIONIZE WELLNESS WITH FORA® 6 Connect

VITALITY

WELLNESS

ESSENTIAL

Discover the pinnacle of health management with FORA® 6 Connect. Experience cutting-edge technology as it accurately measures six vital parameters, ensuring that you stay well-informed. Seamlessly manage your data through Bluetooth connectivity, revolutionizing how you approach your health.



Hassle-free Data
Management

HEALTHY LIVING WITH FORA[®] 6 *Connect*

FORA[®] 6 Connect offers cutting-edge technology to measure six vital health parameters.



Blood Glucose
Hematocrit
Hemoglobin



β -Ketone



Uric Acid



Total Cholesterol



↑ ↓
Telehealth
Cloud



High Accuracy
& Precision



Trends &
Day Average



Bluetooth
Connectivity



Free Download


iFORA HM



WELLNESS MADE EASY

Save Time,
Money & Effort

Reduce Your
Carbon Footprint



FORA® 6 Connect

Understand your health
anywhere, anytime.

Complete all tests within just
5 minutes, at your comfort.



Test Strips	BG-HCT-HB	β-Ketone (KB)	Total Cholesterol (TCH)	Uric Acid (UA)
Blood Sample Requirement	0.5 µL	0.8 µL	3.0 µL	1.0 µL
Reaction Time	5 Seconds	10 Seconds	60 Seconds	15 Seconds

*All FORA® 6 Test Strips Series can be purchased separately.

FORA® 6 Connect



Model No.

FORA® 6 Connect

Benefits

FORA® 6 Connect is a multi-parameter monitoring system allowing testing of **Blood Glucose (BG)**, **Hematocrit (HCT)**, **Hemoglobin (HB)**, **β-Ketone (KB)**, **Uric Acid (UA)**, and **Total Cholesterol (TCH)**.

Measuring Range (BG test)

10 ~ 600 mg/dL (0.55 ~ 33.3 mmol/L)

Meal Tag

AC, PC, AutoQC

LCD Backlight

Yes

Strip Indication Light

Yes

Strip Ejection

Yes

Connectivity (external output)

Bluetooth

Power Source

One 1.5V AAA alkaline battery

Battery Life

1,000 times

Memory Capacity

1,000 measurement results with respective dates and time

Dimensions

H 89.8 x W 54.9 x D 18 mm

Weight

46.1g (without batteries)



- High accuracy test strips made with gold
- Advanced 5-electrode Technology
- Advanced Superior Sip-In Technology



- Sync to the free iFORA App and check your measuring records anytime and anywhere.



- Save Time and Money with the hospital-grade FORA® 6 Connect and enjoy unparalleled efficiency in your health management.

FORA®
ForaCare Europe
www.foracare.de

Contact us: +49 (0)621 76021 410

Instagram: [fora_lifestyle](#)

Instagram: [foracare_de](#)

Facebook: [FORA Deutschland](#)

Importer and Distributor:

Smart OTC GmbH

Markircher Straße 9A,
68229 Mannheim,
Deutschland

TEL: +49 (0)621 76021410

FAX: +49 (0)621 76021444